BEAUTIFLY

B-Derma
PRO
series

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Technical data B—Derma PRO

charging type	wireless charging
power voltage:	5V—1A
battery:	3.7V 1500mAh /1.85Wh
charging duration:	24hours
usage time:	≥ 1.5h
size:	65mm*40mm*40mm
material used:	ABS MS
suction power:	61Kpa±5Kpa
idle time:	90 days (if the unit hasn't been used for more than 90day's, please fully charge it before use)
included:	device, wireless charge base, 7 heads, 6 sponges, charging cabel, user manual

Product introduction

1. Replacement probe
2. Battery display
3. Levels indication
4. Working display
5. Operating time (5mins)
6. Levels button
7. ON/OFF/Pause button
8. Wireless charge

Microcrystalline tip

It can be used on the whole face with the exception of the eyes and in areas on which hard skin such as the knees or elbows were formed.

Big Round Hole Tip

Except eyes. Can be used on any part of the face.

Small Round Hole Tip

Can be used on sensitive skin around the eyes, and also other parts of the face.

Oval Hole Tip

Can be used at eye corner, nose-wing and skins around the mouth.

Working principle

Vacuum suction

The vacuum massage system improves the firmness and oval of the face. Pressure massage promotes better circulation of blood and lymph, making the skin tense, better supplied with blood, and swelling and sagging skin reduced.

Microcrystalline head

It produces over 100,000 natural abrasive particles that make the dead skin thoroughly removed from the face. The skin remains smooth and delicate.

Blue photon light

Intended for acne-prone but sensitive skin. Bio-effects: it has a positive effect on acne skin, reduces inflammation, regulates sebaceous glands and the amount of sebum secreted, the skin becomes more tense and enlarged pores become less visible.

Using area

Microcrystalline head: can be used on any part of the face except close to the eyes. It can also be applied to body areas covered with hard calloused skin, e.g. on elbows or knees.

Large round hole beauty head: can be used on any part of the face except close to the eyes

Small round hole beauty head: can be used on sensitive skin around the eyes, also on other parts of the face.

Oval hole beauty head: can be used on the face, nose and other parts of the face and neck.

Operation instructions

1. Press the power button for 2 seconds to turn on the device, press the power button for 2 seconds during use to turn it off.
2. The white light will have a flashing during charge, once fully charged the light will stop flashing and the battery will light red.
3. To pause the device, hold the power button temporary. The battery will flash if the energy level is low. The device will shut down automatically after 5 minutes of inactivity. Turn on the unit if you need to use it again.
4. When the device is working, the timer will stop when it runs for 5 minutes. If it needs to be used again, please restart it. After the pause, the time will be recalculated.
5. After using a certain level to shut down, reboot will be in the original gear.

Application

- Discoloration
- Scars
- Stretch marks
- Extended, very visible pores
- Comedones and acne (avoid inflammatory changes, phlegias)
- Wrinkles
- Skin devoid of a healthy color with the first signs of photoaging

How to use

Microdermabrasion

1. Please thoroughly wash your face. This will makes your pores open up.

Please do not apply the device at the same spot on your skin more than 3 seconds. It needs to move up and down constantly.

Please attention that sunscre moisturizing in time to avoid skin hurt after using the product, it is

3. recommended to apply toner or cold water on the face, this will close the pores.
4. Please adjust the suction pressure based on your skin type.
5. After used, it is normal that your face to be a little red and sensitive, the redness will disappear in 30 min.
6. Please do not use device in conjunction with other beauty tools or remedies.

Pore cleansing

1. Use blackhead leading-out liquid or apply face with hot towel for 3-5 minutes to open the pores and export the blackheads.
2. Put on the big round hole tip or small round hole tip.

Please point the tip on the blackhead or acne directly to suck, you can choose the suction level according to your skin conditions.
3. After sucking, pat cold water or apply mask to shrink pores.

When sucking the blackhead, you should move slowly on the blackhead area and repeat this action in one direction instead of moving in circles. As for acne, please suck it then release it immediately.

Face fine lines removing

1. Put the oval hole tip on the device.
2. Use the device to massage eye corner, nose-wing, lips. You can choose the suction level according to your need.

You can suck and release skin parts that need caring quickly. Repeat operation steps mentioned above to improve blood circulation, increase skin elasticity and prevent wrinkles.

Lift skin and smooth away wrinkles with low-gracie suction. Please don't suck the same part for too long as that may cause redness on your skin. If in this case, please stop using temporarily.

Contraindications

- Active acne phlegosis
- Rosacea
- Therapy with retinoids
- Acidic, lasers - at least 4 weeks before microdermabrasion
- Psoriasis
- Virotic, bacterial and fungal diseases of the skin vascular skin
- Vascular skin
- Scarring
- Distent continuity of the skin (fresh wounds)
- Skin cancer or suspected lesions in the place undergoing surger
- Fillers in the face - at least 4 weeks before microdermabrasion, otherwise, the preparation may move.

Frequency of treatment

The procedure is performed approximately once every two weeks in a series of approximately 6-12 treatments at intervals of 7 to 28 days, depending on the expected effect. About 3-4 days after the treatment there is a micro exfoliation. It does not require the use of an ampoule or serum for the procedure.

Post-treatment care

Immediately after microdermabrasion the skin may be visibly red and slightly irritated. The mechanical exfoliation treatment of stratum corneum can be combined with iontophoresis or needdeless mesotherapy. Thanks to this, the active ingredients absorability at these treatments is much greater and more beneficial than using, for example, needdeless mesotherapy.

It is necessary to apply creams with a filter immediately after the treatment, preferably du ring the whole series. Avoid chlorinated water, sauna, solarium, exposure to the sun.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out basket on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may lead to health detriment or environmental pollution during further disposal or recycling of the equipment.

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Beautifly Sp. z o.o., Plac Bankowy 2, 00-095 Warszawa

English

Hereby, Beautifly Sp. z o. o. declares that the device equipment Beautifly B-Derma PRO is in compliance with Directive 2014/30/UE.

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