BEAUTIFLY

Derma
Blush
series

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- www.beautifly.eu
- www.youtube.com/beautiflyeu

Tension:	5V
Current:	450mAh
Battery:	Lithium-ion battery
Power:	1.6 W
Size:	163.5x52x38mm
Libra:	205g
Suction power:	50-60 kPa

Security:

Before using the product, read this user manual and follow the recommendations contained therein when using the device.

Device description [page 2-f1]:

1. Charging indicator
2. Suction power indicator
3. Power button/suction power adjustment
4. USB port

Description of heads:

1. Small round suction head. Low suction power For people who have dry or sensitive skin and those who are just starting treatments with the device.
2. Oval suction head. Medium suction power/special shape. For lifting and firming treatments.
3. Round suction head. Gentle/medium suction power. For people who have oily skin, acne or blackheads.
4. Small microcrystalline suction head
5. Oval microcrystalline suction head
6. Round microcrystalline suction head

Service:

Enable/Disable

Press and hold the power button for 3 seconds, the product is turned on, the display is on and the suction power displays 1 level (soft).

Power regulation:

Short press the power button to adjust the suction power. Each short press increases the suction power by one level, in the order: soft > normal > strong > ultra strong.

Application:

- Discoloration
- Scars
- Stretch marks
- Enlarged, very visible pores
- Blackheads and acne (avoid inflammatory and purulent lesions)
- Wrinkles
- Skin devoid of healthy color with the first signs of photoaging

Contraindications:

- Active purulent acne
- Rosacea
- Retinoid therapy
- Acid treatments, lasers - at least 4 weeks interval before microdermabrasion
- Psoriasis
- Viral, bacterial and fungal skin diseases
- Vascular skin
- Tendency to scarring
- Broken skin continuity (fresh wounds)
- Skin cancer or suspicious moles in the area being treated
- Facial fillers - at least 4 weeks before microdermabrasion, otherwise the preparation may move.
- Do not wash the device with water, the device is not waterproof.

Turn the suction head to a 90-degree position and slowly move it over your face. Don't hold in the same part for a long time. (Before removing the suction head from the skin, tilt it to the side to release the pressurized air. This will allow you to remove it safely without irritating the skin. Do not pull or move it directly to avoid pulling the skin and causing redness.)

More information about the company and products [page 2 - f2]:

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[EN] Rules for handling electrical and electronic equipment: Electrical and electronic equipment as well as batteries and accumulators must not be combined with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Dispose of waste electrical and electronic equipment as well as batteries and accumulators at your local separate collection point or to your dealer. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Incorrect handling or damage to them may lead to damage to health or environmental pollution when the equipment is further disposed of or recycled.

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[EN] Summary declaration of conformity: The manufacturer hereby declares under his sole responsibility that the above product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of legislation, Member States in the field of electromagnetic compatibility and Directive 2011/65/EC of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment and Commission Delegated Directive (EU) 2015/863 of 2011. On 31 March 2015, amends Annex II to Directive 2011/65/EU of the European Parliament and of the Council on the list of restricted substances.