BEAUTIFLY

Derma Glow

PRO

series

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- www.beautifly.eu
- www.youtube.com/beautiflyeu

Rated voltage:	3.7V
Rated power:	< 1.6W
Battery capacity:	400mAh
Size:	175mm x 46mm x 55.5mm
Weight:	305g
Working time:	> 100 minutes
Charging time:	2 hours
Set includes:	4 interchangeable heads, USB cable, case, replaceable filters, user manual

Safety

Read this manual before using the product and follow the instructions in it when using the device.

Product description [page 2 - f1]

1. Head
2. Display screen
3. Power button
4. USB charging port
5. Back charging indicator

Replaceable suction heads [page 2 - f2]

1. Small round suction head - Facilitates the removal of blackheads from hard to reach places.
2. Round suction head - Vacuum lifting produces a gentle, stimulating massaging effect that improves microcirculation and tightens the skin.
3. Elliptical suction head - It can be used at the eye corner, nose-wing and skins around the mouth.
4. Microcrystalline round suction head - Removes dead and calloused epidermis, acne scars.
5. Microcrystalline tower suction head
6. Microcrystalline small round suction head

Charging

1. Charge the device for 8 hours at most.
2. Make sure the device is off when charging.
3. You can connect the device to a power adapter via USB cable, or directly to a computer.
4. When charging, the back charge indicator flashes, and when fully charged, the charge indicator lights up.

Precautions

- Dry your hands before using the device.
- Do not rinse or immerse the device under water.
- Do the daily skin care after using this device.
- Do not use other devices for mechanical facial cleansing while using

this device.

- Do not put too much pressure on your skin while using, as the device itself has a strong suction.
- The suction heads should be disinfected to avoid the emergence and spread of viral, bacterial and fungal diseases of the skin.

Contraindications

- Active purulent acne
- Rosacea
- Retinoid therapy
- Acid treatments, lasers - at least 4 weeks before microdermabrasion
- Psoriasis
- Viral, bacterial and fungal diseases of the skin
- Vascular skin
- Tendency to scarring
- Broken skin (fresh wounds)
- Skin tumors or suspicious moles in the treated area
- Fillers in the face - at least 4 weeks before performing microdermabrasion, otherwise the preparation may move.
- Do not wash the device with water, the device is not waterproof.

More information about the company and products [page 2 - f3]

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[EN] Rules for handling electrical and electronic equipment: Electrical and electronic equipment as well as batteries and accumulators must not be combined with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Dispose of waste electrical and electronic equipment as well as batteries and accumulators at your local separate collection point or to your dealer. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Incorrect handling or damage to them may lead to damage to health or environmental pollution when the equipment is further disposed of or recycled.

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[EN] Summary declaration of conformity: The manufacturer hereby declares under its sole responsibility that the above product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances.

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