BEAUTIFLY

Scrub Mini

Young
series

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- www.beautifly.eu
- www.youtube.com/beautiflyeu

Rated voltage:	AC 100V~240V 50HZ/60HZ
Output voltage:	5V/400mA
Battery:	lithium battery
Weight:	70g
Ultrasound:	25KHZ
Power:	2W
Size:	140*42*20mm
Kit Includes:	user's manual, USB cable, accessories

Safety

Please read this manual before using the product and follow the instructions in it when using the device.

Device description [page 2 - f1]

1. Stainless steel blade
2. Protective cover
- 3 Power button/mode switch
4. USB eye charging port

Product features

- Provides deep cleansing of the skin.
- Reduces pigment spots and breaks down pigment.
- Smoothes wrinkles and improves skin elasticity.
- Improves blood circulation and lymphatic metabolism.
- Increases the regenerative capacity and permeability of skin cells.
- Improves the condition and appearance of the skin.
- Is easy to carry and navigate.

Charging

Note: Fully charge the battery before using the device for the first time.

- Check the condition of the lithium battery before charging.
- Connect the USB charging cable to the USB interface.
- Connect to the USB interface of your computer, an adapter (sold separately) or the lithium battery of your mobile phone.
- The LED flashes continuously during charging and lights up when charging is complete.

Notes on charging

- Charge the device if it does not display a start-up status when the power button is pressed for a long time.
- During operation, the device will sound a warning tone and then the shutdown state will start. This indicates that the battery level is insufficient. The product must be recharged immediately.
- After 10 minutes of use, the product automatically enters standby mode and, if not used for one minute, switches off automatically.
- It takes approximately 120 minutes to fully charge. (It takes more than 5 hours to charge before first use).
- Once the battery is fully charged, the device can be used for approximately 60 minutes.
- It is recommended to recharge the device periodically.

Operating instructions

Important: The switch lights flash when changing modes.

Cleaning mode

1. Press the power button for 1.5 seconds to switch the unit on.
2. The first mode - cleansing will start automatically when the unit is switched on.
3. Daily care: reduction of blackheads, whiteheads and acne. The ultrasonic vibrations penetrate deep into the muscles, deeply cleansing the skin. Use together with a facial cleanser.

Moisturising mode

1. Press the power button for 1.5 seconds to switch on the device.
2. Press the power button once to start the second mode - Moisturising.
3. Ultrasonic vibrations and negative ion make the active ingredients in cosmetics penetrate deep into the skin.

Lift mode

1. Press the power button for 1.5 seconds to switch on the device.
2. Press the power button twice to activate the third mode - Lifting.
3. Micro-currents generated by ultrasound and EMS stimulate cell activity. This process stimulates collagen and elastin production, resulting in firmer and smoother skin.

Cleaning and maintenance

- Protect the device from moisture, do not place it in the dishwasher or microwave.
- Switch off the device before cleaning.
- Do not use harsh, caustic or volatile thinner cleaning materials.
- Use a soft cloth or cotton swab moistened with cold or warm water to gently wipe the unit.
- Do not rinse the device under running water.

More information about the company and products [page 2 - f2]

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[EN] Rules for handling electrical and electronic equipment: Electrical and electronic equipment as well as batteries and accumulators must not be combined with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Dispose of waste electrical and electronic equipment as well as batteries and accumulators at your local separate collection point or to your dealer. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Incorrect handling or damage to them may lead to damage to health or environmental pollution when the equipment is further disposed of or recycled. The device was placed on the market after August 13, 2005

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[EN] Summary declaration of conformity: The manufacturer hereby declares under its sole responsibility that the above product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances.

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